

AED Inspection

School:	Date:
AED Make/Model:	Serial Number:
Inspector:	

Step	Tasks	N/A	Pass	Fail
Pads Expiration Date Pediatric Expiration Date	Verify pads are in good condition with no visible damage to cable, packaging, or plug. Verify expiration date. Verify pediatric pad condition and expiration date if applicable.			
Battery Status Expiration Date	Verify battery status and expiration date. Replace secondary 9V battery if applicable and older than 6 months.			
Housing	Verify housing is in good condition with no visible damage to handles, latches, and audible alarm, if applicable.			
Supplies	Verify the needed supplies are present, including scissors, gloves, breathing device, alcohol swab, and means for hair removal.			
General Condition/Visible Damage	Verify that no visible damage is noted on the AED, including body; screen, if applicable; buttons, plugs, etc. (Describe in comments.)			
Function	Verify the unit functions properly: powers on and off properly, gives audible prompts, and has no error codes or other discrepancies.			
Comments				

